

The UV Cube: High-Capacity Disinfection Helping Keep Staff and Patients Safe



Safe for All Environments



The UV Cube KR-31XL-CART



US Pentagon-Selected Technology

Proven High-Capacity UV-C Infection Prevention

The UV Cube—AUVS' largest Germicidal Enclosure—features a 13-cubic-foot chamber to disinfect large items and meet the demands of high-throughput environments.

Disinfect multiple laptops, phones, tablets, keyboards, backpacks, purses, briefcases, hand tools, personal items and PPE including N95 and KN95 masks (during shortages) in just seconds. Simply place items inside, without stacking, lock the door and press a single button to begin the disinfection cycle.

The UV Cube is based on AUVS' patented technology developed for the US Pentagon and now used throughout more than 700 Hospital Groups.

Featuring:

- Proven 254nm UV-C Disinfection
- Fast 90-Second Cycle Time
- One-Button Ease of Use
- Independent Lab-Tested Disinfection:
 - Human Coronavirus - 99.999%
 - MRSA – 99.99%
 - C.Diff Spores 99.99%
 - Norovirus – 99.999%
- Electronics-Safe Technology Proven to Not Damage Your Devices
- Chemical Free Operation, Producing No Waste Products



Disinfect up to
16 Chromebooks
in just **90 seconds.**

Made in the USA / EPA Est No. 99021-GA-1

Help Control Outbreaks and Prevent Infection.



Distributed by...
Hospital Safety Solutions
844-889-6757
www.hospitalsafetysolutions.com

